


# Planmeca ProSensor HD

Intraoral sensor





# Planmeca ProSensor HD

## Unparalleled Image Quality

The advanced imaging sensor with a fibre-optic layer offers real HD image quality.

True resolution of over 20 lp/mm.\*

Theoretical resolution 33 lp/mm.



\* High resolution mode only

# Planmeca ProSensor HD

## Unparalleled Image Quality



Sharp, low-noise and high contrast images.



The wide dynamic ensures successful results every time.

# Planmeca ProSensor HD

## Imaging from Children to Adults

To fulfill all your intraoral imaging needs, the sensor is available in three sizes (S0, S1, S2).

Rounded sensor edges make the procedure comfortable for the patients.



# Planmeca ProSensor HD

## Imaging from Children to Adults

White sensor color enhances visibility to ease positioning.

Images captured in seconds.



**PLANMECA**

# Planmeca ProSensor HD Easy One Hand Use

Magnetic connector is easy to attach using only one hand.

Ensures that the connector is always on correct position.



# Planmeca ProSensor HD Instant Visual Feedback

Control box is equipped with colour-coded LED light providing instant visual feedback of the imaging procedure





# Planmeca ProSensor HD

## Easy Infection Control

Sensors are hermetically sealed and can be fully immersed in disinfectant.





# Planmeca ProSensor HD

## Carefree Choice

Sensor cable has been reinforced and includes only two wires.

The five-year warranty program ensures a safe investment.



# Planmeca ProSensor HD

## Integration

USB version also available

Mac OS support

Seamless integration with  
Planmeca ProX intraoral x-ray  
unit.



# Planmeca ProSensor HD

## Romexis Software

Native **Mac OS X** support and Windows support

**Complete** imaging system for 2D and 3D images

**Versatile modules** for implant planning, Cephalometric analysis, Airways, TMJ analysis, Smile Design and more.

### Compatibility

Full DICOM support

PMbridge, TWAIN, Quick Launch

Free Romexis **Viewer** application

**Planmeca mRomexis** for mobile use.

**Romexis Cloud** service to share images online

